

TEST NAME	TEST CODE	COLLECTION TUBE	VOLUME (mL)		TURN-AROUND APPROX.
			REQUIRED	MINIMUM	
DNA Deletion Duplication Array	DDDA	EDTA (purple top)	4	2	4 to 6 weeks
Quantitative PCR	qPCR	EDTA (purple top)	4	2	4 to 6 weeks

SPECIMEN COLLECTION / SUBMISSION / DOCUMENTATION / REQUISITION FORM

- 1** Specimens must be received by the Advanced Genomics Laboratory (AGEN) within 3 days from time of blood draw.
- 2** Specimens are to be peripheral blood collected in EDTA.
- 3** Proper labeling of specimen tube:
 - Patient's First and Last Name (correctly spelled)
 - Date of Birth
 - Medical Record # (Unique Identifier)
 - Date Drawn
- 4** Identical information must be provided between the specimen collection tube(s) and the submitted requisition form.
- 5** The lab does not bill third parties; therefore, please be sure to provide complete institutional billing information.

SHIPPING AND HANDLING

- 1** Specimens are to be shipped Sunday through Thursday for first priority overnight earliest delivery by 8AM.
- 2** Please fax tracking # and name of delivery service to Advanced Genomics Laboratory at 414-955-6128.
- 3** DOT and IATA requirements for shipping diagnostic proper labeling of specimen tube:
 - Place specimen in a leakproof primary receptacle.
 - Place primary receptacle in leakproof secondary vessel.
 - Place absorbent material between the primary and secondary packaging.
 - An itemized list of contents must be enclosed between the secondary packaging and the outer packaging.
 - Package and Air Waybill must show the text: "Diagnostic Specimen" and labeled with UN3373 designation.
- 5** Specimens must maintain room temperature during transport; recommend insulated shipping containers.
- 6** All delivery carrier (local and non-local), deliver to the Medical College of Wisconsin LOADING DOCK.
 - LOADING DOCK entrance is located off of Watertown Plank Road between 87th & 92nd streets.
 - Medical College of Wisconsin LOADING DOCK hours are 8AM to 3PM (Monday thru Friday); closed on weekends.

SPECIMEN DELIVERY ADDRESS

Send samples at room temperature to:

Medical College of Wisconsin
 Advanced Genomics Laboratory (AGEN)
 Attn: Rachel Lorier
 TBRC / CRI Rm C2388
 8701 Watertown Plank Road
 Milwaukee, WI 53226
 Website: www.mcw.edu/AGEN
 Email: www.AGEN@mcw.edu

Phone: 414-955-2358

Fax: 414-955-6128